
Enabling a Healthier World

Lonza

Moving Towards Sustainable BET Testing

PyroGene® rFC Assay

Katrin Pauls | 22 August 2024

Business Use Only





Speaker:

Katrin Pauls

Market Development and Scientific Affairs Manager

Katrin studied Microbiology and Virology at the University in Duesseldorf, Germany and has over 15 years' experience in QC Microbiology for Pharma, Food and aviation industries. In her current role as Market Development and Scientific Liaison Manager for Testing Solutions at Lonza's Bioscience Cologne site Katrin is responsible for developing markets and leading communicating on technical aspects of sustainable testing solutions.



Traditional Bacterial Endotoxin Tests (BET)

- PYROGENT® Gel Clot Assays
- PYROGENT® 5000 Turbidimetric Assays
- Kinetic-QCL® Chromogenic Assays



Sustainable Pyrogen & Endotoxin Tests

- PyroGene® Recombinant Factor C (rFC) Assays
- PyroCell® Monocyte Activation Test (MAT) Rapid Systems



Instrumentation, Software, & Accessories

- Multimode, Absorbance and Fluorescence Readers
- PyroTec® PRO Automated Robotic Solution
- WinKQCL® Endotoxin Analysis Software
- Accessories and Consumables



Training & Support

- e-Learning Certification Modules
- QC Insider® Toolbox
- Testing Services
- Field Support Services
- Scientific Support

Lonza has been providing endotoxin detection solutions and services since the 1970s



The horseshoe crab – a “sustainable” creature

- > By fossil records over 445 million years old
- > Occupied a niche in the ecosystem, 4 species are still around today
- > Mandatory actions to stabilize the population, e.g.
 - Quota for Fishing/ Bait industry
 - Measures against erosion, habitat shifting and alteration due to coastline development
 - Actions against climate change and severe weather
 - **Use of recombinant methods**



ASMFC¹ – *Limulus polyphemus* (LAL)

- > Overall neutral trend in population.
- > Multiple challenges from pressure on habitats – action plans for protection

IUCN² – *Tachypleus gigas*, *T. tridentatus* and *Carcinoscorpius rotundicauda*

- > Endangered and data deficient
- > No action plan in place

1. <http://www.asmf.org/species/horseshoe-crab>
2. <https://www.iucnredlist.org>


Recombinant Factor C Assay – A Sustainable BET Test?

A simple answer... Yes!




- Eliminates the dependency on LAL and TAL as a natural resource
- Secures the supply chain – addresses the constant growth in the number of annual BET tests
- Absence of biological variability, higher lot-to-lot consistency
- No false positive results - Glucan path eliminated
- Comparability of LAL/TAL and rFC demonstrated by a large number of studies

Endotoxin Detection



Recombinant Factor C Assay - Supported by Regulations




European Pharmacopeia Compendial	United States Pharmacopeia → Compendial (New)	Japanese Pharmacopeia Alternative	Chinese Pharmacopeia Compendial
<ul style="list-style-type: none">> Ph. Eur. 2.6.32 general analytical chapter> Ph. Eur. 5.1.10 Guidance for using BET tests> New (2025): Ph. Eur. 5.1.13 Guideline for selecting a pyrogenicity test (→ monograph methods)	<ul style="list-style-type: none">> USP <1085> general chapter> New (2025): USP <86> general analytical chapter "recombinant methods"	<ul style="list-style-type: none">> JP G4-4-180 BET & alternative methods using recombinant proteins	<ul style="list-style-type: none">> ChP 1143 Guideline for BET applications
 <p><i>Guidance for Industry: Questions and Answers</i> "If a manufacturer chooses to use a recombinant factor C based assay method validation should be in accordance with the requirements of USP <85> , "BET" and USP <1225>, Validation of compendial procedures."</p> <p style="text-align: right;"><small>*Pharmacopeial Forum 36:1 Jan/Feb 2010</small></p>			

History of PyroGene® Recombinant Factor C Assay




PyroGene® rFC Assay history




History of PyroGene® Recombinant Factor C Assay (2)



PyroGene® rFC Assay history



PyroGene® Recombinant Factor C Assay – Global Footprint



Customer Locations:

- | | |
|----------------|--------------------------|
| Australia | Italy |
| Austria | Japan |
| Belarus | Luxembourg |
| Belgium | Mexico |
| Brazil | Netherlands |
| Canada | Norway |
| China | Poland |
| Croatia | Saudi Arabia |
| Czech Republic | Singapore |
| Denmark | South Korea |
| Estonia | Spain |
| France | Switzerland |
| Germany | Taiwan |
| Hong Kong | Thailand |
| Hungary | Turkey |
| India | United Kingdom |
| Ireland | United States of America |

The sustainable Bacterial Endotoxins Test

- > A quantitative, single-step, endpoint photometric method with sensitivity down to 0.005 EU/mL
- > Assay uses liquid reagents, takes 1 hour at 37° C
- > Follows BET methodology and guidelines, the assay is performed in a standard 96-well microplate
- > The assay read-out is performed in a Fluorescence Reader (PyroWave®, Nebula® → Ex/Em 380/440 nm) controlled by WinKQCL® Software



The sustainable Bacterial Endotoxins Test

- > Uses standards calibrated to RSE to produce a standard curve
- > Routine testing is as per quantitative methods
- > Uses a PPC at or near the middle of the standard curve
- > Assessment of assay data as per USP/EP/JP
- > Follows validation protocols as described
- > The primary difference between the rFC and LAL methods is that the rFC reagent is not sourced from HSC blood



PyroGene® Recombinant Factor C Assay - QuickGuide




This is a step-by-step guide depicting how to perform the PyroGene® Recombinant Factor C assay. Prior to initiating the assay procedure, allow reagent vials to equilibrate to room temperature. The incubating microplate reader should also be turned on and a plate template created in the WinKQCL™ Software.

Step 3

Label the tubes with the appropriate endotoxin concentration and add 0.75 ml of LRW to the 5.0 EU/ml tube. Add 0.9 ml of LRW to each of the remaining tubes.

(Example based on a test with an operating standard curve of 0.005–5.0 EU/ml.)



Step 5

Dispense 100 µl of the LRW blank, endotoxin standards, product samples, positive controls, etc. into the appropriate wells of the microplate.



Step 1

Reconstitute Control Standard Endotoxin (CSE) with LAL Reagent Water (LRW) to yield a solution containing 20 EU/ml.



Step 2

Vortex for 15 minutes.



Step 4

Prepare a series of endotoxin standards.



Step 6

Pre-incubate the plate for ≥ 10 minutes at 37°C ± 1°C in the microplate reader.



Step 7

During incubation period, prepare the working reagent, which consists of fluorogenic substrate, assay buffer and rFC enzyme solution in a 5:4:1 ratio, respectively. Add enzyme last to the buffered substrate.



PyroGene® rFC Assay is a standard photometric assay, remarkably like kinetic chromogenic or kinetic turbidimetric LAL



rFC is the same binding protein operating in the LAL assay



The activated recombinant Factor C enzyme cleaves a substrate directly instead of activating another enzyme in a series (the LAL cascade)



The substrate has a fluorescent tag, which gives a wide dynamic range with better resolution

Phased Approach for Adoption



Bacterial Endotoxin Testing (BET) is an essential quality control test for raw materials, in-process samples and final product release



When adopting a recombinant assay into your endotoxin testing workflow a conservative, phased approach could be considered



Phase 1 – raw material and in-process samples (such as in-process water). Other samples such as dialysate and research samples should be considered



Phase 2 – final products



Phase 3 – automation



Raw/In-Process Materials



Final Product




Automation

PyroGene® rFC Assay Validation

Validation Protocol



Please contact your local scientific support team for validation protocol and support!



[Validation Strategy for New Recombinant Factor C Users](#)


Method Comparability – PyroGene® Assay Validation



A Recombinant Factor C Procedure for the Detection of Gram-negative Bacterial Endotoxin

Pharmacopeial Forum, Vol. 36 (1), Jan.-Feb. 2010

- > Multi-center study – sites in US and Europe
- > Set-up followed requirements of USP Chapter <1225>
- > Ten different pharmaceuticals or related products
- > Tested at MVD, 1/2 MVD and 1/10 MVD
- > **The rFC method was comparable to LAL procedures but offered more specificity (non-reactive to glucans)**




Results of Pyrogene® rFC Assay and LAL Assays Are Comparable



Comparison of four endotoxin detection reagents in measuring naturally contaminated endotoxin levels in four representative parenteral products

DocID: GUID-39D7842E-76C1-4CF2-B2F6-A10B2857B74A_10101_en-US

- > Four parenteral products (Saline, Insulin, Acyclovir, Gentamicin) spiked with natural water to simulate a breach in the purification system
- > Products were tested using a kinetic turbidimetric method, a kinetic chromogenic method, and two rFC methods
- > All three assay types safely detected the contaminants
- > **There is no “best method” for every parenteral product**



Why PyroGene® Recombinant Factor C Assay?



In Summary

- > Lonza pioneered the commercialization of an equivalent, reliable recombinant endotoxin detection method
- > The PyroGene® rFC Assay has been gaining acceptance world-wide
- > With the recombinant Factor C assay, you can meet the demands of supply security while protecting a natural resource
- > A phased implementation of the rFC assay offers a conservative approach
- > Reliable lot-to-lot consistency, statistically more robust than LAL methods
- > Adopting recombinant reagents provides many benefits and further coupling those reagents with automation helps to reduce errors in the BET assay process and prepare for the future

PyroGene® rFC Assay – Kit Size



Catalogue Number	No of tests/reactions	Components	Sensitivity (EU/mL)
50-658U	192 tests	2 x 96 tests/vial rFC Enzyme 2 x 6mL vials Fluorogenic Substrate 2 x 5mL vials rFC Assay Buffer 2 vials E. coli Endotoxin 2 x 30mL vials LAL Reagent Water	0.005 to 5
50-658NV	2880 tests	30 x 96 tests/vial rFC Enzyme 30 x 6mL vials Fluorogenic Substrate 30 x 5mL vials rFC Assay Buffer 10 vials E. coli Endotoxin	0.005 to 5






Become a QC Insider

QC Insider Toolbox - Endotoxin Expertise at Your Fingertips®



An online suite of comprehensive resources from beginner to advanced support tools, training modules, videos, webinars and a library of information that helps QC professionals perform classic and recombinant bacterial endotoxins test.

Visit: www.lonza.com/qcinsider

<p>Support →</p> <p>Access to guidance and detailed information about software support, recertification and more.</p> <ul style="list-style-type: none">• Recertification Services• Reader Install & Maintenance• Software/IT Support• Scientific Support• Testing Services• Automation 	<p>Training →</p> <p>Self-directed resources that will help you build on your bacterial endotoxins test expertise.</p> <ul style="list-style-type: none">• e-Learning Modules• Videos• Webinars• Seminars 	<p>Library →</p> <p>Technical resources that will help ensure your success with the bacterial endotoxins test.</p> <ul style="list-style-type: none">• Certificates of Analysis• Package Inserts• Quick Guides• Scientific Posters• Tech Tips• White Papers• Assay Selection Guide• Presentations 
---	---	---

Questions?



Any questions?

Please submit your questions using the questions pane at the bottom of your screen

Interested to learn more?



Visit our website

bioscience.lonza.com



Follow us on LinkedIn®

www.linkedin.com/showcase/lonza-qc-testing-solutions



Subscribe to our eNews

www.lonza.com/enews

All trademarks belong to Lonza, registered in USA, EU or CH or to third party owners and used only for informational purposes. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products.

For more details: www.lonza.com/legal. ©2023 Lonza. RT-PP189 02/23

Lonza

Thank you!

All trademarks belong to Lonza, registered in USA, EU or CH or to third party owners and used only for informational purposes. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. For more details: www.lonza.com/legal
© 2024 Lonza

Katrin Pauls | 22 August 2024

