


# Sustainable Pharmaceutical Quality Control with Recombinant Factor C Technology

Ganesh Chaudhari





# AGENDA

- Bacterial Endotoxin Testing – Challenges & Implications
- Introduction of Recombinant Factor C (rFC)
- The ENDONEXT™ Range – How to do thing differently?
- Summary and Q&A

# AN HISTORIC FAMILY COMMITMENT TO MEDICINE AND PUBLIC HEALTH WORLDWIDE

- **A Pasteurian tradition:** Marcel Mérieux worked with Louis Pasteur in 1894.
- **A commitment that transcends four generations:** ever since the creation of Institut Mérieux, a generation has worked after another to expand its legacy.



**Marcel Mérieux**  
1894 - Student of **Louis Pasteur**



**Marcel Mérieux**  
1897 - Creation of Institut Mérieux




**Dr Charles Mérieux**  
1937 - Dr Charles Mérieux took up the reins



**Dr Alain Mérieux**  
1963 - Creation of bioMérieux

**Alexandre Mérieux**  
2017 - Chairman & CEO of bioMérieux

# INSTITUT MÉRIEUX: SERVING PUBLIC HEALTH FOR OVER A CENTURY



\* Companies which are no longer part of the Institut Mérieux since 1994 and where created by the Mérieux family

# DEVELOPMENT OF ENDOTOXIN TESTING

## Limulus Amebocyte Lysate (LAL) Test

- Introduced in 1970s
- Detects endotoxins from GNB (in vitro)
- Principle: gelation, colour or turbidity change
- Requires bleeding protected horseshoe crabs



## Monocyte Activation Test (MAT)

- Introduced in 2010 in Ph. Eur. 2.6.30
- Non-animal based in vitro pyrogen test
- Represents a full replacement of the rabbit test
- ELISA assay to measure cytokine release from human monocytes in response to pyrogens



## Rabbit Pyrogen Test (RPT)

- Introduced in 1942, the original pyrogen test
- Detects endotoxin and non-endotoxin pyrogens (in vivo)
- Principle: rise in body temperature of rabbits following IV injection
- Not sustainable and has limitations



## Recombinant Factor C (rFC) Test

- First registered in 2003 (PyroGene)
- Specific to endotoxin from GNB
- Principle: fluorescence signal
- Non-animal based & sustainable alternative to LAL/TAL

## What's next?



# ENDOTOXIN TESTING CHALLENGES



## Cost of Reagents & Equipment

Does the method fit in my budget? How do I do the testing without compromising product and patient safety?



## Sustainability compliance

Does my method require animal based raw materials? Does it help with sustainability initiatives of my business? Are my supply chain practices sustainable?



## Regulatory / Validation

Is it recognised by regulatory authorities? Is it compendial or alternative? Who does the PQ? Can I justify the spend?



## Vendor Support

Are they experts, can I rely on them? Do they have sufficient stock levels?



## Sample Throughput

How many samples can I run? How much time does it take? Is the method automated or manual?



## Sample Matrix Interference

Does the method work for my samples? Do I need to use dilutions with water or buffers?



## Training / Skills Required

How easy is it to train my operator? Do I need additional resources? What are the chances for errors?




# INTRODUCTION OF RECOMBINANT FACTOR C



# FACTOR C IS THE SPECIFIC SENSOR

## LAL Cascade




LAL is a biological enzyme system, and no 2 lots will react the same



BIOMÉRIEUX

## Recombinant Factor C (rFC)




Exact copy of Factor C produced synthetically

No extraction from the Limulus blood

Simple Reaction without cascade

No  $\beta$ -glucan / sugars interference



rFC shows superior batch consistency





# rFC MEETS PHARMACOPEIA CONTROLS FOR BET

Solution	Endotoxin concentration	Solution to which endotoxin is added	Number of replicates
A	None	Test solution	Not less than 2
B	Middle concentration of the standard curve	Test solution	Not less than 2
C	At least 3 concentrations (lowest concentration is designated $\lambda$ )	Water for BET	Each concentration not less than 2
D	None	Water for BET	Not less than 2




max 21 samples


# TESTING WITH RFC IS JUST LIKE LAL

No surprise – same known requirements

Sensitivity:	Assay time:
0.05 EU/mL	20 min
0.005 EU/mL	60 min
0.001 EU/mL	120 min



Fluorescence



**EU/mL**  
Endotoxin units per millilitre



bioMérieux - ENDONEXT  
Report - Assay 2/23/2022 9:34:36 AM  
Reader name: 11212121  
Assay Status: VALID  
Assay template: E23\_Method Validation\_0209  
Client version: 3.11.13  
Reader mode: E23E17X  
Assay started at: 2/23/2022 9:34:14 AM  
Page 1/12

**Final Report**  
Assay 2/23/2022 9:34:36 AM

**Summary results**

Sample name	ID	Lot Number	Release Limit	Result	Status
WEF	1	12345	N/A	<0.005 EU/ml	VALID
WEF	2	12345	N/A	<0.005 EU/ml	VALID
WEF	3	12345	N/A	<0.005 EU/ml	VALID
WEF	4	12345	N/A	<0.005 EU/ml	VALID
WEF	5	12345	N/A	<0.005 EU/ml	VALID
WEF	6	12345	N/A	<0.005 EU/ml	VALID
WEF	7	12345	N/A	<0.005 EU/ml	VALID
WEF	8	12345	N/A	<0.005 EU/ml	VALID
WEF	9	12345	N/A	<0.005 EU/ml	VALID
WEF	10	12345	N/A	<0.005 EU/ml	VALID
WEF	11	12345	N/A	<0.005 EU/ml	VALID
WEF	12	12345	N/A	<0.005 EU/ml	VALID
WEF	13	12345	N/A	<0.005 EU/ml	VALID
WEF	14	12345	N/A	<0.005 EU/ml	VALID
WEF	15	12345	N/A	<0.005 EU/ml	VALID
WEF	16	12345	N/A	<0.005 EU/ml	VALID
WEF	17	12345	N/A	<0.005 EU/ml	VALID
WEF	18	12345	N/A	<0.005 EU/ml	VALID
WEF	19	12345	N/A	<0.005 EU/ml	VALID
WEF	20	12345	N/A	<0.005 EU/ml	VALID
WEF	21	12345	N/A	<0.005 EU/ml	VALID

# THE ENDONEXT™ RANGE

## HOW TO DO THING DIFFERENTLY?



# THE ENDONEXT™ PORTFOLIO


## REAGENTS

Wide applications; from rapid, high-throughput water testing (ENDOZYME II GO) to complex samples (ENDOLISA) & Low Endotoxin Recovery (ENDO-RS)




## INSTRUMENT

BioTek Synergy HTX multimode reader (fluorescence & absorbance) with integrated incubation and shaking function, and a dual reagent injector module



## SOFTWARE

ENDONEXT™ Software is an integrated, 21 CFR Part 11 compliant solution for data analysis, data management & full reporting



## AUTOMATION

Semi- and fully-automated platforms for QC laboratories with higher throughput - in collaboration with INTEGRA & TECAN



## ACCESSORIES



ENDOGRADE® certified Water & Glass Test Tubes



# EZIGO : FLEXIBILITY, SAVINGS & PERFORMANCE

Manual processing >>

Semi-Automation >>

Full-Automation

## Standard Microplate ENDOZYME II GO


Standard Reconstitution

Prep. of Standard Curve  
+ vortex + deposit

Preparation of samples

Adding PPCs

Filling of Plate



Reagent preparation

Addition of reagent

Successive Reads and Results


Not needed

Not needed

Preparation of samples

Not needed

Filling of Plate




Reagent preparation

Addition of reagent

Successive Reads and Results


## INTEGRA Productive Pipetting



## TECAN



Hands-on time



Time to result  
(1 reader)

1 x  = 40 min

1 x  = 22 min

6 x  = 60 min

# RESULTS IN HALF THE TIME



## ENDOZYME® II GO – The rapid endotoxin test



CES Standards  
50 - 0.005 EU/mL

Negative controls  
(Blank)

PPC Controls

	1	2	3	4	5	6	7	8	9	10	11	12
A	STD 50	STD 50	PPC SPL 1	PPC SPL 2	PPC SPL 3	PPC SPL 4	PPC SPL 5	PPC SPL 6	PPC SPL 7	PPC SPL 8	PPC SPL 9	PPC SPL 10
B	STD 5	STD 5	PPC SPL 1	PPC SPL 2	PPC SPL 3	PPC SPL 4	PPC SPL 5	PPC SPL 6	PPC SPL 7	PPC SPL 8	PPC SPL 9	PPC SPL 10
C	STD 0.5	STD 0.5	SPL 1	SPL 2	SPL 3	SPL 4	SPL 5	SPL 6	SPL 7	SPL 8	SPL 9	SPL 10
D	STD 0.05	STD 0.05	SPL 1	SPL 2	SPL 3	SPL 4	SPL 5	SPL 6	SPL 7	SPL 8	SPL 9	SPL 10
E	STD 0.005	STD 0.005	SPL 11	SPL 12	SPL 33	SPL 14	SPL 15	SPL 16	SPL 17	SPL 18	SPL 19	SPL 20
F	BLK	BLK	SPL 11	SPL 12	SPL 33	SPL 14	SPL 15	SPL 16	SPL 17	SPL 18	SPL 19	SPL 20
G	BLK	BLK	PPC SPL 11	PPC SPL 12	PPC SPL 13	PPC SPL 14	PPC SPL 15	PPC SPL 16	PPC SPL 17	PPC SPL 18	PPC SPL 19	PPC SPL 20
H	PPC Control	PPC Control	PPC SPL 11	PPC SPL 12	PPC SPL 13	PPC SPL 14	PPC SPL 15	PPC SPL 16	PPC SPL 17	PPC SPL 18	PPC SPL 19	PPC SPL 20

PPC  
0.5 EU/mL

Samples  
(empty wells)

PPC  
0.5 EU/mL



# ENDOLISA® – A SOLUTION FOR COMPLEX MATRICES

ENDOLISA® revolutionizes endotoxin testing of complex samples **with a unique built-in sample preparation**



Broad measurement range:  
**0.05 – 500 EU/ml**

- **ELISA-like format** featuring a 96-well plate **pre-coated** with a specific **endotoxin-binding phage protein**
- **Overcomes limits** of traditional methods such as inhibition and enhancement
- Intended for **complex samples** with interfering substances present in the product matrix
- Important assay for **Endotoxin Demasking** in combination with ENDO-RS® (Low Endotoxin Recovery / Hold Time Studies)
- Described in the **European Pharmacopeia** as solution for removing interferences (Chapter 5.1.10, section 9)



# [ REGULATORY STATUS OF RFC WORLDWIDE ]



# REGULATORY LANDSCAPE OF rFC



## US Pharmacopoeia

<86> Bacterial Endotoxins Test Using Recombinant Reagents

**Compendial**  
PROPOSED



## EU Pharmacopoeia

Ph. Eur. 2.6.32 Test for Bacterial Endotoxins with Recombinant Factor C (rFC)

**Compendial** 2021

Inclusion of rFC in the monographs:

- Water, purified (0008)
- Water for injections (0169)
- Gene Therapy & mRNA Products (proposed)



Ministry of Food and Drug Safety

## Korean Pharmacopoeia

General Information Test No. 23 Endotoxin Testing Method Using Recombinant Factor C

**Alternative** 2023



## Japanese Pharmacopoeia

G4-4-180 Bacterial Endotoxins Test & Alternative Methods using Recombinant Protein-reagents for Endotoxin Assay

**Alternative** 2021



## FDA

Guidance for Industry, Pyrogen and Endotoxins Testing: Q&A

MAPP 5310.7

Supports and accepts the use of rFC to replace LAL



British Pharmacopoeia

## British Pharmacopoeia

Added rFC to Appendix XIV C. Test for Bacterial Endotoxins

**Compendial** 2024



## Indian Pharmacopoeia

Guidelines on the BET Introduction of Alternate Test Methods

**Alternative**



## Chinese Pharmacopoeia

G4-4-180 Bacterial Endotoxins ChP 9251 Guideline for BET Application (introduction of rFC)

**Alternative** 2020



## EAEU Pharmacopoeia

To add a chapter on Test for Bacterial Endotoxins with Recombinant Factor C (rFC) based on Ph. Eur. 2.6.32

**Compendial**  
PROPOSED

# SUMMARY & KEY TAKEAWAYS




# EMBRACING RFC FOR LAB TRANSFORMATION



## EFFICIENCY

### *Time and Cost Savings*

- Decrease in retests & OOS investigations
- Standardized results (No variability due to animal variability)
- Reduce false positives (Not sensitive to Factor G)
- Consistent results batch-to-batch
- Easy staff training: more time for added value tasks in the lab



## TRACEABILITY

### *Standardized Results for Data Integrity and fewer risks of error*

- Reduced hands-on time
- Reduced manual pipetting-related errors
- Increased sample throughput
- Improved assay performance by better reproducibility & accuracy
- Tracking and trending of endotoxin results in a compliant software



## SUSTAINABILITY

### *HSC-Free for Ethics & Compliance*

- Reduction of waste & flexibility
- Meet 3R principles: animal free
- Supply chain sustainability:
  - no production risk linked to ecological matters
  - easy to ship around the world





PIONEERING DIAGNOSTICS

BIOMÉRIEUX